



YOUR WIDEX HEARING JOURNEY

JUNE 2026

WIDEX

SOUND LIKE NO OTHER

YOUR LISTENING AND COMMUNICATION NEEDS



AI-enhanced conversations with **Clarity Boost** - Exclusively available in Widex Allure AI RIC

	TV AND PHONE STREAMING	CONVERSATION IN QUIET	COMFORT DURING TRANSPORT	LISTENING TO MUSIC	COMFORT OUTDOORS	CONVERSATION DURING TRANSPORT	CONVERSATION IN MODERATE NOISE	CONVERSATIONS OUTDOORS	CONVERSATION IN LOUD NOISE (RESTAURANT OR PARTY)
	Stay connected to speech or music streamed from your TV or smartphone.	A balance of clear speech and accurate environmental sounds for a natural sound experience.	Keep noise comfortable when traveling by car or public transportation.	Enjoy music with a natural and authentic sound quality.	Hear comfortably outdoors, even in windy conditions.	Enhance conversations when traveling by car or public transportation.	Hear speech above a moderately noisy background, like a family dinner.	Hear speech while outdoors in windy conditions.	Hear speech in noisy, crowded environments, such as restaurants or a party.
440	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
330	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	
220	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★			
110	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★				
50	★★★★★	★★★★★	★★★★★	★★★★★					

YOUR PERSONALISED HEARING JOURNEY

1

Understand your hearing loss with the help of a qualified Independent Hearing Care Professional and discuss **your audiogram** results.

2

Explore **your listening & communication needs** based on your individualised lifestyle and choose your technology performance level.

3

Choose your style of hearing aid that is best suited to you.


4

Choose your accessories and connectivity options to maximise your communication experience.

YOUR AUDIOGRAM

Client's Name

Clinician's Name



An audiogram is a graph which gives a detailed description of your hearing abilities for different frequencies. The lower the curve, the louder the sound must be in order for you to hear clearly. Looking at this graph, your hearing care professional can tell whether you are suffering from hearing loss, and if so, how severe it is.

CHOOSE YOUR STYLE

RECHARGEABLE

RECEIVER-IN-CANAL (RIC) MODELS



AI RIC R D
• AI Clarity Boost
• Smartphone Connectivity



RIC R D
• Smartphone Connectivity



SmartRIC
• Smart Design
• Smartphone Connectivity



BTE R D
• Smartphone Connectivity



RIC 312 D
• Smartphone Connectivity
• Size 312 Battery



RIC 10
• Smallest RIC
• Size 10 Battery



BTE 13 D
• Smartphone Connectivity
• Size 13 Battery



BTE 312
• Slim BTE
• Size 312 Battery

RECHARGEABLE CUSTOM MODEL



ITE R D
• Smart Phone Connectivity

ITE COLOURS



Black



Mocha



Espresso



Beige



XP/ITC/ITE
• Size 312 Battery
• Telecoil (optional)
• Wireless Programming



CIC/CIC-M
• Size 10 Battery
• CIC Wireless Programming
• CIC-M Non-wireless Programming

COLOURS

Choose between 14 different colors for RIC/BTE models.



Mediterranean
Turquoise



Deep
Blue



Shocking
Pink



Sporty
Red



Titanium
Grey



Chestnut
Brown



Golden
Brown



Autumn
Beige



Pearl
White



Tech
Black



Silver
Grey



Dark
Cherry



Honey
Blonde



Rose
Gold

	Mediterranean Turquoise	Deep Blue	Shocking Pink	Sporty Red	Titanium Grey	Chestnut Brown	Golden Brown	Autumn Beige	Pearl White	Tech Black	Silver Grey	Dark Cherry	Honey Blonde	Rose Gold
ALLURE RIC/BTE					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MOMENT RIC/BTE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
MOMENT SmartRIC										✓	✓	✓	✓	✓

CHOOSE YOUR ACCESSORIES

TV SOLUTIONS



TV PLAY 2 TELEVISION STREAMER

Streams high quality sound from your TV wirelessly to Bluetooth hearing aids.



TV DEX

Streams TV to hearing aids wirelessly. Analogue only (digital-to-analogue converter required for modern TVs).



COMMUNICATION AND STREAMING SOLUTIONS



SOUND ASSIST

A 5-in-1 multi functional device that enhances communication in complex listening situations. Bluetooth hands-free communication, Remote control & telecoil function.



SOUNDCONNECT

Direct streaming stereo high-quality audio from personal computers/laptops.



COM-DEX

Hands-free device that streams high quality audio from most Bluetooth devices to your wireless hearing aids.



RC-DEX

Compact user-friendly mini remote control for discrete volume/program changes.



CHARGERS



RIC/BTE STANDARD CHARGER

Compact and easy to use standard charger.



RIC/BTE CHARGE N CLEAN

Take better care with your premium 3-in-1 charger.



RIC PORTABLE CHARGER

With a portable charger, you will enjoy hassle-free charging, mobility, and peace of mind.



Qi wireless charging



ITE CHARGER

This rechargeable device brings all the comfort of a custom device, without fiddling with tiny batteries.



SMARTPHONE APP OPTIONS



WIDEX ALLURE™

- Control device functions including directionality
- Access AI Quick Assistant
- Compatible with all Bluetooth hearing aids



WIDEX MOMENT™

- Control device functions including directionality
- Access MySound*
- Compatible with all Bluetooth hearing aids

*available in 220 and above



WIDEX ZEN TINNITUS

- Tools & information to help manage tinnitus



WIDEX TONELINK

- Discreetly change programs, volume & directional focus
- Available with all Widex models
- Compatible with all Bluetooth hearing aids



WIDEX MAGNIFY™

- Control device functions including directionality
- Compatible with all Bluetooth hearing aids



Scan the QR code to check your phone compatibility

 .5ms

The faster an amplified sound gets to your ear, the more accurately it is presented to the brain.

Widex delivers sound fastest, by far¹

Why does that matter?

 95%

of wearers agree the sound is natural and clear with greater satisfaction participating in daily life²

How can that help?

 10x

as many wearers are very satisfied in noisy backgrounds using PureSound™ with ZeroDelay™ compared to other hearing aids³

What else can Widex offer?

 80%

of wearers with tinnitus experienced a reduction of tinnitus severity using Widex SoundRelax™⁴

DANISH DESIGNED & ENGINEERED

At Widex, we believe that in most situations natural sound is better than aggressive processing. However, in specific noisy situations where you need support, we use technology to reduce background noise, helping you focus on what you want to hear.



¹Slugocki C, Kuk F, Korhonen P, Ruperto N. Neural encoding of the stimulus envelope facilitated by Widex ZeroDelay technology. Hearing Review 2020;27(8):28-31.

²Balling LW, Townend O, Helmink D. Sound quality in real life-Not just for experts. Hearing Review 2021;28:27-30.

³Balling LW, Townend O, Helmink D. Sound quality for all: The benefits of ultra-fast signal processing in hearing aids. Hearing Review 2021;28(9):32-35.

⁴Balling LW, Caporali S, Parker D, Jeppesen AM, Helmink D. Expanding options for effective tinnitus management. Hearing Review 2023;30(2):12-16.